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Humidification Water Set

Water Container Set

for CO₂ Incubators with Active Humidification

1. Scope of delivery

1.1 Tubular bag set

There are two optional bag sets with all necessary connecting parts:

- Art. no. 8009-1083 Set with 3 empty tubular bags:
- Art. no. 8009-1058 Set with 12 empty tubular bags (annual demand)

Each bag set includes:

- 3 or 12 empty tubular bags (art.no. 6002-0626. 6002-0626)
- 3 two-piece tube connectors with 3-5mm nozzles (art.no. 6009-0322)
- 3 two-piece tube connectors with 5-7mm nozzles (art.no. 6009-0323)
- 3 or 12 gas fine filters (sterile filter) (art.no. 6014-0012)
- 6 m silicone tube (DN 3mm) (art.no. 6008-0267)
- 1m filling tube (DN 10 mm) (art.no. 6008-0268)
- Instructions (art.no. 7001-0385)

1.2 Humidification water set

Optional humidification water set with all necessary connecting parts:

- Art.no. 8009-1104 Set with 10 one-liter bottles, prefilled with sterile, deionized water (annual demand)

The humidification water set includes:


- 10 one-liter bottles, prefilled with sterile, deionized water (art.no. 1007-0314)
- 3 spike connectors (transfer sets) for sterile liquids (art.no. 6007-0131)
- 3 two-piece tube connectors with 3-5mm nozzles (art.no. 6009-0322)
- 2 m silicone tube (DN 3mm) (art.no. 6008-0267)
- 3 gas fine filters (sterile filter) (art.no. 6014-0012)
- 20 cm silicone tube (DN 6 mm) (art.no. 6008-0009)
- Instructions (art.no. 7001-0385)

By using the humidification water set you will achieve unsurpassed contamination safety of your samples in the BINDER CB CO₂ incubator with active humidification.

- (1) Pump head
- (1a) Pump head cover
- (1b) Tube retainer

- (1c) Rotor
- (2) Pump connector

2.2 Water supply via freshwater can

- (12) Pump head
- (12a) Pump head cover
- (12b) Tube retainer
- (12c) Rotor
- (13) Pump connector
- (24) Silicone tube DN 3mm to water can

3. Refilling the freshwater bag with the filling hose

	<p>For easier handling, you can temporarily release the pump head from the pump connector on the incubator.</p>
---	---

The Bag Set contains a filling tube (DN 10 mm) of 1 m length.

- Remove the connection between both parts of the tube connector
- Attach the filling tube to the tube connector part which is still on the freshwater bag
- Fill the freshwater bag. Observe water quality and temperature.
- Remove the filling tube from the tube connector on the freshwater bag
- Reassemble the two parts of the tube connector.


4. Replacing the silicone tube of the fresh water supply system

6 m silicone tube DN 3mm are included with the Bag Set. To operate the freshwater supply a total length of 32 cm is required.

- Turn the pump head (1) counterclockwise until it stops and remove it from the pump connector (2).
- Pull up the tongue of the tube retainer (1b) up to “Click” and remove it. Take the rotor (1c) out of the pump head cover (1a). Remove the silicone tube.
- Lead the new silicone tube along the rotor. Insert the rotor with the tube into the pump head cover. Keep the hose tight and do not pinch it.

The end of the tube for connecting to the incubator should be about 6 cm long. With a total length of the silicone tube of 32 cm, the tube end for connecting the freshwater bag has a length of approx. 18 cm

- Press the tube retainer (1b) into the pump head cover until it snaps into place

Observe the correct orientation of the connections!


- Insert the pump head (1) into the pump connector (2) and turn it clockwise until it snaps into place.
- Attach the silicone tube to the incubator connection
- Place the freshwater bag with the opening to the left in the recess behind the chamber door.

5. Connection / change of the tubular freshwater bag

5.1 Replacing the tube connectors for freshwater supply

3 two-piece tube connectors with 3-5mm nozzles und 3 two-piece tube connectors with 5-7mm nozzles are included with the Bag Set.





State of delivery:

Two-piece tube connector with 3-5mm nozzles
(Art.no. 6009-0322)

Two-piece tube connector with 5-7mm nozzles
(Art.no. 6009-0323)

For connection to the silicone tube and the freshwater bag, the hose connectors must be combined accordingly:

To silicone tube: 3-5 mm nozzle		To freshwater bag: 5-7 mm nozzle	
Combination 1:	Opening NW 7 mm		 Opening NW 4 mm
Combination 2: (recommended)	Opening NW 4 mm		 Opening NW 7 mm

5.2 Connecting the freshwater bag

Put the filled freshwater bag on the nozzle of the hose connector.


6. Connection / change of the prefilled water bottle

The humidification water set includes all necessary connecting parts.

Place the tube connectors on the silicone tube of the water system.

Remove the protective cap of the spike connector (transfer set) and connect the nozzle of the tube connector with the spike connectors using a part of the DN 6 mm silicone tube.

Then remove the second protective cap of the spike connector and insert the piercing pin into the water bottle through the silicone seal.


7. Connecting / refilling the freshwater can

7.1 Connecting the tubing

Mount the freshwater can: Attach the magnetic can holder to the left chamber side and insert the freshwater can.

Then you can establish the connection between pump and freshwater can.

(Image shows single chamber connection)

- | | |
|--|---|
| (17) Freshwater can | (22) Reducing connector |
| (18) Magnetic can holder (water can support) | (23) Wire tension clamp |
| (19) Hose coupling with aqua stop | (24) Silicone tube DN 3mm to the pump (or via T-type tube connector (T-piece) to another chamber) |
| (19a) Removal button | (12) Pump |
| (20) Angular coupling plug with aqua stop | |
| (21) Silicone tube d 10 | |

Water supply for a single chamber:

- Plug the silicone tube DN 3mm (24) coming from the pump (12) on the reducing connector (22)

Water supply for several chambers:

You can branch off the water supply via T-pieces for additional chambers.

- Divide the silicone tube DN 3mm (24) into pieces of a suitable length. Take into account the height of the magnetic can holder
- Put a piece of the silicone tube DN 3mm (24) on the reducing connector (22) and connect a T-piece to the free end
- If there are more than two chambers, you can connect another T-piece using a piece of the silicone tube DN 3mm (24)
- For each chamber, plug the silicone tube DN 3mm (24) coming from the pump (12) onto the free outlet of the T-piece

7.2 Refilling the freshwater can

There are two possibilities:

- Remove the freshwater can from the water can support. Press the grey removal button on the hose coupling, then the angle coupling plug can be removed. A small amount of water can escape (filling amount of the angular coupling plug). The water can outlet remains closed due to the aqua stop.
- Remove the lid and fill the can on site with a suitable aid (water hose, watering can). When filling on site, the existing hose connection can remain on the freshwater can.

8. Replacing the gas fine filter (sterile filter) for ambient air

The air tube for inlet of the ambient air is located in the upper right corner in the recess for the freshwater bag.

We recommend replacing the gas fine filter (sterile filter) on this tube about once a month.

To do this, hold the sterile filter and pull out a bit the air tube connected to it. Unscrew the filter and screw on a new filter. Then push the tube back in.

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meets the provisions of the directive 93/42/EEC in consideration of the change of the directive 2007/47/EEC and its transpositions in national laws which apply to it.

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**TÜV Rheinland LGA Products GmbH