

Stacking options for BINDER CO2 incubators, C and CB series

CHAMBER WITH ROLLER BASE	C 150/170			CB (E3/E6)			CB (E7)		
Chamber model	C 150 (E1/E2)	C 170 (E3)	C 170 (E7)	CB 53/60	CB 150/160	CB 210/220	CB 56	CB 170	CB 260
Height of chamber [mm]	820	870	870	720	920	1070	670	870	1020
Height of roller base [mm]	118	118	118	118	118	118	118	118	118
Height of chamber on roller base [mm]	938	988	988	838	1038	1188	788	988	1138
Art. No. Roller base	9051-0024	9051-0029	9051-0028	9051-0031	9051-0028	9051-0029	9051-0043	9051-0028	9051-0044

Roller base



Roller base with model CB 170



STACK WITH STACKING ADAPTER	C 150/170			CB (E3/E6)			CB (E7)		
Upper chamber model	C 150 (E1/E2)	C 170 (E3)	C 170 (E7)	CB 53/60	CB 150/160	CB 210/220	CB 56	CB 170	CB 260
Lower chamber model	C 150 (E1/E2)	C 170 (E3)	C 170 (E7)	CB 53/60	CB 150/160	CB 210/220	CB 56	CB 170	CB 260
Height of chamber [mm]	820	870	870	720	920	1070	670	870	1020
Height of stacking adapter [mm]	89	89	89	89	89	89	89	89	89
Height of two chambers with stacking adapter (mm)	1729	1829	1829	1529	1929	2229	1429	1829	2129
Art. No. Stacking adapter	9051-0025	9051-0036	9051-0035	9051-0037	9051-0035	9051-0036	9051-0038	9051-0035	9051-0039


Stacking adapter



Stacking adapter with two C 170



STACK WITH ROLLER BASE AND STACKING ADAPTER	C 150/170			CB (E3/E6)			CB (E7)		
Height of roller base [mm]	118	118	118	118	118	118	118	118	118
Height of two chambers with stacking adapter on roller base [mm]	1847	1947	1947	1647	2047	2347	1547	1947	2247
Art. No. Roller base	9051-0024	9051-0029	9051-0028	9051-0031	9051-0028	9051-0029	9051-0043	9051-0028	9051-0044



Roller base and Stacking adapter with two C 170



STACK ON STACKING FRAME WITH CASTERS	C 150/170			CB (E3/E6)			CB (E7)		
Upper chamber model	-	-	-	-	CB 150/160	CB 210/220	-	CB 170	CB 260
Lower chamber model	-	-	-	-	CB 150/160	CB 210/220	-	CB 170	CB 260
Height of two chambers on stacking frame [mm]	-	-	-	-	2054	2354	-	1955	2255
Art. No. Stacking frame	-	-	-	-	9051-0020	9051-0023	-	9051-0040	9051-0041

Stack on stacking frame with casters



Stack on stacking frame with casters with model CB 170 und C 170



Further stack combinations can be requested by contacting BINDER INDIVIDUAL: [go2binder.com/en-BINDER-INDIVIDUAL](https://www.binder.com/en-BINDER-INDIVIDUAL)