








Water supply system

		Series KBF	Series KBF P	Series KBF LQC	Series KMF	Series KBF-S/ KBF-S ECO	
Standard:	Water supply via pressurized, deionized water External	•	•	•	•	–	
Accessory:	BINDER PURE AQUA SERVICE System for preparation or complete desalination of tap water, complete set containing PURE AQUA 300 single-use cartridge, measuring device, and all necessary connecting parts	•	•	•	•	•	
Accessory:	Water supply set Flexible water and waste water management. External, for hanging from the back of the device	•	•	•	•	–	
Accessory:	Water supply set Flexible water and waste water management. External, free-standing	•	•	•	–	–	
Standard:	Convenient water supply Independent water supply via tank	–	–	–	–	•	 
Accessory:	Flexible positioning of the water container with magnetic holder, hose adapter for water container with practical quick coupling, and wa- ter stop function. External, with magnetic holder	–	–	–	–	•	

• available | – not available